

Dotmatrix Printer for Weighing Balance DP

SALIENT FEATURES

HSN Code - 90160020



- ◆ Dot Matrix printing Method
- ◆ Paper width - 58mm
- ◆ Baud Rate (bps): 1200, 2400, 4800, 9600, 19200
- ◆ Buttons and Indicator Lights: provide LF and SEL buttons and indicator light (green)
- ◆ Interface RS232C
- ◆ Power Supply - DC 5V(±5%),3A
- ◆ Date & Time - NA, If the connected instrument has date and time functionality, it can print the date and time along with the other data.



DP Dotmatrix Printer for Weighing balance



Dp5- Cable for printer - M2535M
DP6 - Cable for printer - M3535M
DP7 - Cable for printer - F3535M
DP8 - Cable for printer - F2535M

DP5

TAB4T, TAB9T, TAB61T, KAB1L, KAB2L, KAB4L, KAB7L, LMB16, LMB17T LPB25T, LPB26T, PGB3010, PGB301, PGB200, PFB60, PFB100, PFB300, TTB3, TTB15, TTB31, HPB3000, HPB3010, PGB630, MAB282, MAB250, MAB220LCD, HPB60H, PGB505, PGB205

DP6

TAB2T, TAB3T, KAB5L, KAB6L, EPB6T, EPB7T, EPB8T, PGB1210, LPB23, LPB24, MAB301, MAB201, HPB201, HPB501, HPB1001

DP7

TAB1T, KAB10L, LMB1, LMB11R, LMB20, EMB11, PGB220, PGB321, PGB3110, PGB6100, PGB600, PGB1000, PGB610, TTB20HH, PGB5010M, TTB30M, PGB211M, PGB3001, PGB201D, MAB220T, PGB1MB, PGB640

DP8

KAB8L, KAB9L, TAB10T

STANDARD ACCESSORIES : Main Instrument, Power Cable with Adaptor, RS232C Cable, Printer Paper Roll, Instruction Manual

Data Cable for Weighing Balance to Printer

Dotmatrix Printer for Weighing Balance

DP

HSN Code - 90160020



| Specifications | DP-11, 12, 13, 14 | DP - 15, 16, 17, 18 |
|------------------------------|-------------------|---------------------|
| 128 X 64 Display | ✓ | ✓ |
| Date & Time (built-in clock) | ✓ | ✓ |
| Direct Mode | ✓ | ✓ |
| Weighing Mode | ✓ | ✓ |
| Statistical Mode | - | ✓ |
| Sum Mode | - | ✓ |
| Mean Mode | - | ✓ |
| Formulation Mode | - | ✓ |
| Sample Name Input | - | ✓ |
| Sample Batch number Input | - | ✓ |
| Batch number auto-increase | - | ✓ |
| Data export (use U-Disk) | - | - |
| Data curve print | - | - |
| Barcode Scanner / keyboard | - | - |



E-Catalogue

DOTMATRIX PRINTER PART NO. & SUITABLE MODELS

DP11 - TAB4T, TAB9T, TAB61T, KAB1L, KAB2L, KAB4L, KAB7L, LMB16, LMB17T, LPB25T, LPB26T, PGB3010, PGB200, PGB301, PGB630, PFB60, PFB100, PFB300, TTB3, TTB15, TTB31, HPB3000, HPB3010, MAB282, MAB250, MAB220LCD, HPB60H, PGB505, PGB205

DP12 - TAB2T, TAB3T, KAB5L, KAB6L, EPB6T, EPB7T, EPB8T, PGB1210, LPB23, LPB24, MAB301, MAB201, HPB201, HPB501, HPB1001

DP13 - TAB1T, KAB10L, LMB1, LMB1IR, LMB20, EMB11, PGB220, PGB321, PGB3110, PGB6100, PGB600, PGB1000, PGB610, TTB20HH, PGB5010M, TTB30M, PGB211M, PGB3001, PGB201D, MAB220T, PGB1MB, PGB640

DP14 - KAB8L, KAB9L, TAB10T

DP15 - TAB4T, TAB9T, TAB61T, KAB1L, KAB2L, KAB4L, KAB7L, LMB16, LMB17T, LPB25T, LPB26T, PGB3010, PGB20, PGB301, PGB630, PFB60, PFB100, PFB300, TTB3, TTB15, TTB31, HPB3000, HPB3010, MAB282, MAB250, MAB220LCD, HPB60H, PGB505, PGB205

DP16 - TAB2T, TAB3T, KAB5L, KAB6L, EPB6T, EPB7T, EPB8T, PGB1210, LPB23, LPB24, MAB301, MAB201, HPB201, HPB501, HPB1001

DP17 - TAB1T, KAB10L, LMB1, LMB1IR, LMB20, EMB11, PGB220, PGB321, PGB3110, PGB6100, PGB600, PGB1000, PGB610, TTB20HH, PGB5010M, TTB30M, PGB211M, PGB3001, PGB201D, MAB220T, PGB1MB, PGB640

DP18 - KAB8L, KAB9L, TAB10T

STANDARD ACCESSORIES : Main Instrument, Adaptor, RS232C Cable, Printer Paper Roll, Instruction Manual