



Density Kit for Solid Test
Suitable for WENSAR make Weighing Balance



SDK11, 12, 13, 14, 15, 16
SOLID DENSITY KIT

Suitable for WENSAR make Weighing Balance

A density kit cannot work independently It must be attached with a Analytical / Weighing Balance


Density Kit Pan Holder design perfect fit according Weighing Balance Holder.

HSN : 90160020

Part No. ↓	Density Kit Suitable Model Details
SDK11	TAB1T, KAB4L
SDK12	KAB5L, KAB6L, MAB201, MAB301, HPB201, HPB501, HPB1001
SDK13	TAB2T, TAB3T, EPB6T, EPB7T, EPB8T
SDK14	PGB220, PGB321, PGB211M, PGB201D
SDK15	PGB200, PGB301
SDK16	KAB2L, TAB5T

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument , Mounting Bracket, SS Stand, Glass Beaker, Dropper, Thermometer, Tweezer, Glass Sinker, Hanging Spiral Pan (For weight in Water), Test Certificate.



WDK200, 201, 202, 203, 204, 205
LIQUID DENSITY KIT
 Suitable for WENSAR make Weighing Balance

Density Kit for Liquid Test.
 Suitable for WENSAR make Weighing Balance



A density kit cannot work independently
 It must be attached with a Analytical /
 Weighing Balance

LIQUID DENSITY	
NAME	DENSITY
Water (25°C)	0.997
Water (4°C)	1.000
Alcohol	0.78
Kerosene	0.82
Mercury	13.60


Density Kit Pan Holder design perfect fit according Weighing Balance Holder.

HSN : 90160020

Part No. ↓	Density Kit Suitable Model Details
WDK200	TAB1T, KAB4L
WDK201	KAB5L, KAB6L, MAB201, MAB301, HPB201, HPB501, HPB1001
WDK202	TAB2T, TAB3T, EPB6T, EPB7T, EPB8T,
WDK203	PGB220, PGB321, PGB211M, PGB201D
WDK204	PGB200, PGB301
WDK205	KAB2L

STANDARD ACCESSORIES : Main Instrument , Mounting Bracket, SS Stand, Glass Beaker, Dropper, Thermometer, Tweezer, Glass Sinkers, Test Certificate.

* Specifications are subject to change without Notice



Density Kit for Liquid Test.

Suitable for Other make Weighing Balance



WDK250
LIQUID DENSITY KIT

Suitable for Other make Weighing Balance

WDK250 Density Kit Suitable for other make Weighing Balance, If weighing pan less than 125mm.

If Weighing Balance not having Density function facility, Customer have to calculate manually.

A density kit cannot work independently.It must be attached with a Analytical / Weighing Balance

Density Kit for Liquid Test.

HSN : 90160020

LIQUID DENSITY	
NAME	DENSITY
Water (25°C)	0.997
Water (4°C)	1.000
Alcohol	0.78
Kerosene	0.82
Mercury	13.60

- 1.Hanging basket:
- 2.Density Bracket
- 3.Crystal Ball
- 4.Glass Bowl
- 5.Tank Bracket
- 6.Temperature meter
- 7.Tweezers

STANDARD ACCESSORIES : Main Instrument , Mounting Bracket, SS Stand, Glass Beaker, Dropper, Thermometer, Tweezer, Glass Sinker, Test Certificate.

* Specifications are subject to change without Notice