




SALIENT FEATURES

- ✦ Internal Calibration
- ✦ E.M.F.C Based
- ✦ 7" Touch Screen
- ✦ Halogen Heating System
- ✦ Temperature Range : 40°C ~ 230°C
- ✦ Number Storage up to 200
- ✦ The mode for testing Automatic, Timing, Manual
- ✦ Company Name, Sample Name, Contact No. can print
- ✦ Dual Motor Drive Automatic Open / close
- ✦ Administrator / Operator Password access Login
- ✦ GMP / GLP Compliance


110g / 0.1mg
Moisture Analyzer
EMB11



MENU



PARAMETER SETTINGS



HISTORY RECORD



HISTORY RECORD

Moisture Analyzer

HSN Code - 90160020

SPECIFICATIONS

MODEL	EMB11
Make	WENSAR
Capacity	110g
Readability	0.0001g (0.1mg)
Repeatability	0.0001g
Linearity (±)	0.0003g
Display	7" inches Colour Touch Screen
Display Results	Moisture / Sample Weight / Dry%
Aluminium Sample PAN	Available
Pan Size	100mm
Sensor	EMFC
Calibration	Internal Calibration
Operating Temperature	5°C-35°C
Heating Source	Halogen Lamp
Interface & Communication	RS232C
Downtime Method	Automatic / Manual / Timing
Moisture resolution	0.001%
Temp. Setting	40°C-230°C
Temp. Step	1°C
Time Setting	1-99min
Test Pattern	Standard / Fast / Soft
Storage	Upto 200
Auto Chamber Open/Close	Yes
User Password Protection	Yes
Response Time	≤4Sec
Date & Time	Available
Power Supply	AC 220V / 50-60 Hz
Packing Dimension (LxWxH)	485 x 345 x 355 mm
Gross Weight	10kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Aluminium PAN (50nos), Dust Plate, PAN Handle, PAN Supporter, Power Cable, Instruction Manual, Test Certificate.