





SALIENT FEATURES



- ✦ Capacity - 210g/1mg
- ✦ Loadcell Sensor
- ✦ 7 inches Display
- ✦ Halogen Heating System
- ✦ Temperature Range up to 200°C
- ✦ RS232C Interface
- ✦ The mode for testing Automatic
- ✦ Timing, Manual
- ✦ Wireless Indicator with Rechargeable Battery backup
- ✦ Auto chamber open/ close Function
- ✦ GLP Compliance

210g / 1mg
Moisture Analyzer
LMB20



Moisture Analyzer

HSN Code - 90160010

SPECIFICATIONS

MODEL	LMB20
Make	WENSAR
Capacity	210g
Readability	0.001g (1mg)
Repeatability	0.001g
Linearity (±)	0.003g
Display	7 inches Touch Screen
Display Results	Moisture / Sample weight / Dry%
Aluminium Sample PAN	Available
Pan Size	100mm
Sensor	Loadcell
Calibration	Automatic External
Operating Temperature	5°C - 35°C
Heating Source	Halogen Lamp
Interface & Communication	RS232C & Wireless Indication
Moisture Range	0.00%-100.00%
Moisture Resolution	0.01%
Temp. Setting	40°C~200°C
Temp. Step	1°C
Time Setting	1-99min.
Storage	20 groups
Date & Time	Available
Power Supply	AC 220V / 50-60 Hz
Packing Dimension (LxWxH)	580 x 375 x 455mm
Gross Weight	10kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Indicator, Aluminium PAN, Dust Plate, PAN Handle, PAN Supporter, Power Cable, Instruction Manual, Test Certificate.