



5000g/0.01g
Precision Balance
PGB5010M

SALIENT FEATURES

- ✦ Internal Calibration
- ✦ Touch Screen display (7" full Color HD)
- ✦ Under Weighing
- ✦ Date & Time Setting
- ✦ Piece Counting, Accumulation
- ✦ Check Weighing
- ✦ Percentage Weighing
- ✦ Density Function
- ✦ Net gross weight / sum
- ✦ Overload Alarm
- ✦ RS232C,USB



HSN Code -90160010

SPECIFICATIONS

MODEL	PGB5010M
Make	WENSAR
Capacity	5000g
Readability	0.01g (10mg)
Repeatability	0.01g
Linearity (±)	0.03g
Pan Size	305 x 230mm
Calibration	Internal Calibration
Display	Touch Screen
Operating Temperature	15°C to 35°C
Unit	mg, g, ct, oz, kg, ib, gsm
Application Mode	Piece Counting, Check Weighing, Percentage, Density
Interface	RS232C, USB
Power supply	DC Adaptor, Using AC 220V / 50Hz
Packing Dimension (L x W x H)	460 x 370 x 210mm
Gross Weight	7kg

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Weighing Pan, Power Adaptor, Instruction Manual, Test Certificate.