



SALIENT FEATURES

- ✦ LCD Display with White Back Light.
- ✦ Multifunction Weighing Units.
- ✦ Overload Alarm.
- ✦ In built rechargeable battery.
- ✦ Glass Windshield.
- ✦ RS232C Interface (Optional).
- ✦ Extra Display (Optional).

200g,300g/1mg
Precision Balance
LPB21, LPB22



HSN Code - 90160010

SPECIFICATIONS

MODEL	LPB21	LPB22
Capacity	200g	300g
Make	WENSAR	
Readability	0.001g	
Repeatability	0.001g	
Linearity (±)	0.002g	
Display	LCD Display with white Back light	
Pan Size	80mm Ø	
Tare Range	Full	
RS232C Interface	Optional	
Operating Temperature	10°C to 40°C	
Power Supply	DC Adaptor, using AC 220V / 50-60Hz	
Packing Dimension (LxWxH)	425 x 325 x 425mm	
Gross Weight	5kg	

* Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Adaptor, PAN, Instruction Manual, Test Certificate.