



51g/0.01mg - 220g/0.1mg

100g/0.01mg - 120g/0.1mg

Semi Micro / Analytical Balance

KAB8L, KAB9L

Antivibration Pad - Standard Supply



E - Catalogue

SALIENT FEATURES

- ✦ Electro Magnetic Force Compensation
- ✦ White Backlight LCD Display
- ✦ Detachable Display
- ✦ Fast Sensitive Setting
- ✦ Time Display
- ✦ Weighing Result can Export to Excel
- ✦ Dynamic Weighing Function
- ✦ Pipettes Calibration Function
- ✦ RS232C & USB Interface
- ✦ Piece Counting, Totalizing (Add) Function
- ✦ Automatic Internal Calibration
- ✦ Percentage, Underhook Weighing
- ✦ Multiple Weighing Units

No.	Display	Describe
1	Max: 8888g d= 888.8mg	Model display
2		Internal calibration sign
3	-	Negative weight value
4	0	Stable value
5		Weighing sign
6		Density determination
7		Piece counting
8		Dynamic weighing
9		Summation
10	-----	Progress bar
11	88:88:88	Time setting
12	W1 W2	W1: weight in air, W2: weight in liquid
13		Unauthenticated numbers
14		Unit display



Semi Micro / Analytical Balance

HSN Code - 90160020

SPECIFICATIONS

MODEL	KAB8L	KAB9L
Max. (Capacity)	51g / 220g	100g / 120g
Readability	0.01mg / 0.1mg	
Repeatability	0.01mg / 0.1mg	
Linearity (±)	0.04mg/0.4mg	
Make	WENSAR	
Verification Scale Intervals (e)	1mg	
Calibration	Internal Calibration	
Class	Class I	
Display	White Backlight LCD Display	
Response Time	4 - 6 Sec	
Pan Size (mm)	90 mm Ø	
Tare Range	Full	
Operating Temperature	15°C to 30°C	
Sensitivity Drift (±)	2ppm / °C	
Density	Function Available	
Data Output Mode Selection	Command & Continuous Output	
Key Sound	Yes	
Weighing Chamber Size	162 x 180 x 210mm	
Unit	g, mg oz, ct,	
Application Mode	PCS,%, Weighing, Animal Weighing, Pipettes Calibration, Density, Σ weighing	
Baud rate Options	1200, 2400, 4800, 9600	
Interface	RS232C	
Power Supply	DC9V, Adaptor using AC 220V / 50-60Hz	
Packing Dimension (LxWxH)	490 x 390 x 485mm	
Gross Weight	11kg	

* Specifications are subject to change without notice.

* Density Determination Kit for Liquid or Solid Samples separately (Optional)

STANDARD ACCESSORIES : Main Instruments, Antivibration Pad, PAN, Display Indicator, Power Adaptor, Instruction Manual, Test Certificate.