



SALIENT FEATURES

- ◆ Digital control with LED display.
- ◆ 5 recycle digital timer, 90, 180, 280, 380, 480 seconds.
- ◆ Overheating protection function.
- ◆ Degas function available to improve cleaning efficiency.
- ◆ Built-In cool fan for steady and safe working.
- ◆ Detachable power cord.

Ultrasonic Cleaner

LMUC2D



E - Catalogue

HSN : 84562000

SPECIFICATIONS

MODEL	LMUC2D
Tank Capacity	2.0L
Ultrasonic Power	65W
Heating Power	100W
Degas Function	Yes
Timer	5 cycles: 90s-180s-280s-380s-480s
Ultrasonic Frequency	40 kHz
Tank Material	Stainless Steel SS304
Tank Size (L×W×H)	245×155×80 mm
Overall Size (L×W×H)	305×235×215 mm
Display	Digital control with LED display
Heating Temperature	RT+5°C~ 60°C
Power Supply	AC 220V / 50-60Hz
Packing Dimensions (L×W×H)	345 x 270 x 250mm
Gross Weight	3kg

* Specifications are subject to change without Notice

Note : Ultrasonic Cleaner is not recommended for the following gemstones Opal, Pearl, Emerald, Tanzanite, Malachite, Turquoise, Lapis and Coral.

Standard Accessories : Main Instrument, Power Cord, Basket, Instruction Manual, Test Certificate.

Ultrasonic Cleaners available from 2ltrs to 100ltrs capacity