



Paddle Blender Homogenizer PBH80T

SALIENT FEATURES

- ◆ UV Sterilisation Function
- ◆ Heating Function upto 80°C
- ◆ 4.3" Touch Screen Display
- ◆ 8 Sets of Program can store. adjustable speed & time
- ◆ Glass transparent window is easy to observe & clean
- ◆ Stainless steel with anticorrosive plastic spraying
- ◆ No contact between sample and homogenizer
- ◆ Automatic stop function
- ◆ Front Door & Paddle SS304

Paddle Blender Homogenizer are used mainly in Microbiology application to extract & wash intact Microbes into solution. It is also used in Food Testing Laboratories, Environmental & life Science Research Labs etc..



E - Catalogue

HSN Code : 84798200

SPECIFICATIONS

MODEL	PBH80T
Make	LABMAN
Capacity	4ml ~ 400ml
Temperature range	RT +5°C ~80°C
Timing Setting	Yes
Disinfection wavelength	253.7nm
Slapping speed	3.0~12.0 times of squeezing/second
Speed programmed	8 sets
Adjustable flapping distance	0~50mm
Display	4.3" Touch Screen
Size of sterile homogeneous bag:	170 × 300 mm
Alarm Function, Pause Function	Yes
Paddle Size	145 x 85mm
Power supply	AC 220V / 50Hz
Dimension (LxWxH) & Gross Weight	530 x 330 x 465mm & 25kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Power Cord, Sample Cover 50 nos., Instruction Manual, Test Certificate,