



Ice Flaker - Touchscreen
LMIF30T,50T,100T



Touch operation icon:

- Start-Stop
- Set
- Add
- Reduce
- Timing

SALIENT FEATURES

- ◆ Touchscreen Display
- ◆ Fully Automatic microprocessor control
- ◆ Noiseless operation
- ◆ Compact Design and low maintenance
- ◆ Temperature
- ◆ Front & Back Ventilation design
- ◆ Continuous Ice Flakes Output
- ◆ CFC free compressor
- ◆ Stainless Steel construction
- ◆ Low water level Indication
- ◆ Timer Function
- ◆ SS304



E - Catalogue

HSN Code : 84186990

MODEL	LMIF30T	LMIF50T	LMIF100T
Make	LABMAN		
Capacity of Ice flake (kg per 24 hrs.)	30 kg	50 kg	100 kg
Storage Capacity (kg ~ app.)	12 kg	16 kg	35 kg
Storage bin	ABS Plastic		
Cooling Method	Air		
Ambient temperature	Below 25°C for optimum working		
Construction	Stainless Steel		
Noise level	≤52 dB		
Water Supply	3 / 4 inch BSP		
Refrigerant	R290A		
Input Power	180W	280W	420W
Power supply	AC 220V,50/60Hz		
Dimension LxWxH (mm)	610 x 470 x 780mm	610 x 470 x 820mm	690 x 550 x 900mm
Gross Weight	51kg	52kg	74kg

SPECIFICATIONS

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, Inlet Tube, Drainage Tube, Ice Shovel (Scoop), Instruction Manual, Test Certificate