



**Ice Maker & Crusher**  
**LMIM20C**

## SALIENT FEATURES

- ◆ Ice Crushing Facility.
- ◆ Housing Material Stainless Steel Type.
- ◆ R290A Refrigerant.
- ◆ Water Inlet Manual & Automatic.
- ◆ Bullet Shape Ice.
- ◆ Easy to move through hand grip.
- ◆ Low Water Alarm.
- ◆ Ideal Working Temperature 20°C.
- ◆ Transparent Window on Top Cover to view inside.



E - Catalogue

HSN : 84186990

## SPECIFICATIONS

MODEL	LMIM20C
Ice capacity	20kg/24H
Function	Ice making / Crushing ice
Storage Bin	3kg
Refrigerant	R290A
Ice Shape	Bullet Shape
Housing material	Stainless Steel type
Water inlet method	Manual / Automatic Water inlet
Ice Output Each Cycle	18 pieces
Ice Making Time	10-12 minutes in Ideal Conditions
Input power	180W
Power supply	AC 220V/50Hz
Packing Dimension (L x W x H)	455 x 430 x 460mm
Gross Weight	17Kg

\* Specifications are subject to change without Notice

**Standard Accessories :** Main Instrument, Ice Scoop, Crushed Ice Collecting Cup, Water Inlet Hose, Water Drain Hose, Anti Drain Plug, Inlet Tap Connector, Drain Pipe Connector 1 & 2, Silicon Plug, Instruction Manual, Test Certificate.