



**SALIENT FEATURES**

- ✦ 5L Volume Capacity, PID temperature control Technology
- ✦ LED Display, Stepless Speed regulation
- ✦ Low Noise, Stable Operation
- ✦ PT100 Temperature Sensor
- ✦ High Temperature Resistance
- ✦ Fast Heating Speed
- ✦ Uniform temperature and High efficiency
- ✦ Strong Magnet to ensure Sufficient Suction
- ✦ Heating Method - Electric & Tube Heating

**Oil Bath Magnetic Stirrer**

**LOBMS5**

HSN Code : 84798200



E - Catalogue

**SPECIFICATIONS**

MODEL	LOBMS5
Make	LABMAN
Maximum Volume	5000ml
Display	LED
Speed range	Up to-2800 RPM (Adjustable)
Stirring Power	12 W
Liner Height	145mm
Liner Material	SS304
Temperature Sensor, Control	PT100, PID Temperature
Time Setting	Yes
Max. Temp,	300°C
Temp. Resolution	±1°C
Power Wattage	1000 W
Power supply	AC 220V / 50Hz
Dimension (LxWxH) & Weight	350 x 350 x 340mm & 7kg

\* Specifications are subject to change without notice.

**STANDARD ACCESSORIES** : Main Instrument, Pot Lid 5L, Stir Bar, Temperature Sensor, Temperature Sensor Bracket, Clamp, Beaker Holder, Cross Clamp, Instruction Manual, Test Certificate.