



SALIENT FEATURES

- ✦ Constant Temperature Incubator
- ✦ High precision temperature control, over temperature alarm, inner chamber surround heating, uniform and stable temperature.
- ✦ Compact Design, Light weight, Small volume energy equipped with handle design, Easy to carry & clean
- ✦ Digital Display

Mini Thermostatic Incubator

LMIO6

Applications: It is used for industrial and small enterprises, health, Biochemistry, agricultural science and other Scientific Research Departments to do constant temperature Culture experiments such as bacteria Culture, Breeding, Fermentation etc.



E - Catalogue

HSN Code : 84198990

SPECIFICATIONS

MODEL	LMIO6
Make	LABMAN
Volume	6L
Power	100W
Temp.range	RT +5 ~ 45°C
Ambient temperature	5-30°C
Temp. Resolution & Accuracy	0.1°C & 1°C
Shape size	100 x 125mm
Structure	Portable
Time Setting	Yes
Inner chamber size (LxWxH)	215 x 220 x 205 mm
Outer size mm (LxWxH)	300 x 320 x 355 mm
Power supply	AC 220V / 50Hz
Packing Dimension (LxWxH)	390 x 375 x 415mm
Gross weight	10kg

★ Specifications are subject to change without notice.

STANDARD ACCESSORIES : Main Instrument, SS Tray, Power Adaptor, Instruction Manual, Test Certificate.